

PRACTICAL NO: 3

AIM: Assisting a doctor in diagnosing disease from medical images of patient data

-----CODE-----

```
# !kaggle datasets download -d paulthothymoonney/chest-xray-pneumonia
# !unzip chest-xray-pneumonia.zip

import tensorflow as tf
from tensorflow.keras.preprocessing.image import ImageDataGenerator
import matplotlib.pyplot as plt
import os
import numpy as np
from tensorflow.keras.preprocessing import image

train_dir = "/content/chest_xray/train"
val_dir = "/content/chest_xray/val"
test_dir = "/content/chest_xray/test"

train_datagen = ImageDataGenerator(
    rescale=1./255,
    shear_range=0.2,
    zoom_range=0.2,
    horizontal_flip=True
)

val_datagen = ImageDataGenerator(rescale=1./255)
test_datagen = ImageDataGenerator(rescale=1./255)

train_generator = train_datagen.flow_from_directory(
    train_dir,
    target_size=(150, 150),
    batch_size=32,
    class_mode='binary'
)

val_generator = val_datagen.flow_from_directory(
    val_dir,
    target_size=(150, 150),
    batch_size=32,
    class_mode='binary'
)

test_generator = test_datagen.flow_from_directory(
    test_dir,
    target_size=(150, 150),
    batch_size=32,
    class_mode='binary'
)

model = tf.keras.models.Sequential([
    tf.keras.layers.Conv2D(32, (3, 3), activation='relu', input_shape=(150, 150, 3)),
    tf.keras.layers.MaxPooling2D(2, 2),
    tf.keras.layers.Conv2D(64, (3, 3), activation='relu'),
    tf.keras.layers.MaxPooling2D(2, 2),
    tf.keras.layers.Conv2D(128, (3, 3), activation='relu'),
    tf.keras.layers.MaxPooling2D(2, 2),
    tf.keras.layers.Conv2D(128, (3, 3), activation='relu'),
    tf.keras.layers.MaxPooling2D(2, 2),
    tf.keras.layers.Flatten(),
    tf.keras.layers.Dense(512, activation='relu'),
    tf.keras.layers.Dense(1, activation='sigmoid')
])

model.compile(
    loss='binary_crossentropy',
    optimizer=tf.keras.optimizers.Adam(),
    metrics=['accuracy']
)
```

```

history = model.fit(
    train_generator,
    steps_per_epoch=train_generator.samples // 32,
    epochs=10,
    validation_data=val_generator,
    validation_steps=val_generator.samples // 32
)

test_loss, test_acc = model.evaluate(test_generator, steps=test_generator.samples // 32)
print(f"Test Accuracy: {test_acc * 100:.2f}%")

acc = history.history['accuracy']
val_acc = history.history['val_accuracy']
loss = history.history['loss']
val_loss = history.history['val_loss']

epochs_range = range(len(acc))

plt.figure(figsize=(8, 8))
plt.subplot(1, 2, 1)
plt.plot(epochs_range, acc, label='Training Accuracy')
plt.plot(epochs_range, val_acc, label='Validation Accuracy')
plt.legend(loc='lower right')
plt.title('Training and Validation Accuracy')


plt.subplot(1, 2, 2)
plt.plot(epochs_range, loss, label='Training Loss')
plt.plot(epochs_range, val_loss, label='Validation Loss')
plt.legend(loc='upper right')
plt.title('Training and Validation Loss')
plt.show()

sample_image = r'/content/chest_xray/test/PNEUMONIA/person100_bacteria_475.jpeg'
img = image.load_img(sample_image, target_size=(150, 150))
img_array = image.img_to_array(img) / 255.0
img_array = np.expand_dims(img_array, axis=0)

prediction = model.predict(img_array)
print(f'Predicted Label: {"Pneumonia" if prediction > 0.5 else "Normal"}')

```

-----OUTPUT-----

1/1  1s 1s/step
Predicted Label: Pneumonia